# Bond™ DAB Enhancer

**Catalog No: AR9432**

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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## **Instructions for Use**

Please read before using this product.

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### **Check the integrity of the packaging before use.**

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## Catalog No: AR9432

### Intended Use

*This reagent is for in vitro diagnostic use.*

Bond DAB Enhancer is a ready to use solution that intensifies the DAB staining of immunohistochemical stained, formalin-fixed, paraffin-embedded tissue on the Bond automated system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation).

3,3'-diaminobenzidine (DAB) is a chromogen used for immunohistochemical staining. In the presence of peroxidase enzyme, DAB produces a brown precipitate at the epitope site that is insoluble in alcohol. The Bond DAB Enhancer precipitates as a copper salt at the reaction site producing a darker brown staining deposit.

### Reagents Provided

Bond DAB Enhancer contains 0.5% (w/v) copper sulfate in saline with surfactant.

Total volume = 30 mL.

### Dilution and Mixing

Bond DAB Enhancer is ready to use. Dilution of this reagent is not recommended.

### Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond system.

### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label. Return to 2–8 °C when not in use.

There are no obvious signs indicating contamination and/or instability of Bond DAB Enhancer.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

### Precautions

- This product is intended for *in vitro* diagnostic use.

**BOND DAB Enhancer**                      H412: Harmful to aquatic life with long lasting effects.      P273: Avoid release to the environment.  
Contains Copper Sulphate (<1%).

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.

### Instructions for Use

For use of Bond DAB Enhancer refer to "Dilution and Mixing".

### Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

### Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

### Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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- Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
- Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
- Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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