Bond™ Enzyme Pretreatment Kit

Catalog No: AR9551

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

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Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Catalog No: AR9551

Intended Use

This kit is for in vitro diagnostic use.

Bond Enzyme Pretreatment Kit consists of an enzyme concentrate and an enzyme diluent, which require mixing. The diluted enzyme solution is for enzymatic digestion of formalin-fixed, paraffin-embedded tissue on the Bond automated system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation).

The use of an enzyme pretreatment on formalin-fixed, paraffin-embedded tissue on the Bond automated system exposes epitopes that have been masked by formalin fixation, allowing accessibility of the primary antibody to the epitope.

Reagents and Materials Provided

1. Bond Enzymatic Concentrate (1 mL) contains a proteolytic enzyme (17 mg/mL) and stabilizer.
2. Bond Enzyme Diluent (200 mL) contains Tris-buffered saline, surfactant, and 0.35% ProClin® 950.
3. Three Bond Open Containers (7 mL)

Dilution and Mixing

Dilute before use. The enzyme can be diluted to three different strengths as Enzyme 1, Enzyme 2 or Enzyme 3 in the Bond Open Containers supplied. To make 7 mL of enzyme:

Enzyme 1 - mix 1 drop of Bond Enzyme Concentrate with 7 mL of Bond Enzyme Diluent.

Enzyme 2 - mix 2 drops of Bond Enzyme Concentrate with 7 mL of Bond Enzyme Diluent.

Enzyme 3 - mix 3 drops of Bond Enzyme Concentrate with 7 mL of Bond Enzyme Diluent.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond system.

Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the bottle label.

Diluted Bond Enzyme should be stored at 2–8 °C when not in use, and can be used for the life of the Bond Enzyme Pretreatment Kit.

The signs indicating contamination and/or instability of diluted enzyme solution are turbidity of the solution, odor development, and presence of precipitate.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin® 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- Proteolytic enzyme is harmful. It is irritating to eyes, respiratory system, and skin and may cause sensitization by inhalation. Target organ: Lungs.

BOND ENZYME CONCENTRATE

Contains Proteinase K
(<10%).

GHS08: Health hazard.

Signal words: Danger.

H334: May cause allergy or asthma symptoms or breathing difficulties if inhaled.

P261: Avoid breathing dust/fumes/gas/mist/vapours/spray.

P285: In case of inadequate ventilation wear respiratory protection.

P304+311: IF INHALED: If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing.

P342+311: If experiencing respiratory symptoms: Call a POISON CENTRE or doctor.

P501: Dispose of contents/container to hazardous or special waste collection point.

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.

Instructions for Use

For use of Bond Enzyme Pretreatment Kit refer to "Dilution and Mixing".

Product Specific Limitations

The appropriate enzyme concentration and incubation may vary due to variation in tissue fixation and must be determined by the user. Overdigestion of tissue sections may result in loss of tissue morphology. Negative reagent controls should be used when optimizing retrieval conditions.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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