Bond™ Epitope Retrieval Solution 1

Catalog No: AR9961

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

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Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Catalog No: AR9961

Intended Use

This reagent is for *in vitro* diagnostic use.

Bond Epitope Retrieval Solution 1 is a ready to use epitope retrieval solution for the heat-induced epitope retrieval (HIER) of formalin-fixed, paraffin-embedded tissue on the Bond automated system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation).

The use of a HIER pretreatment on formalin-fixed, paraffin-embedded tissue on the Bond automated system restores epitopes that have been modified by formalin fixation, allowing accessibility of the primary antibody to the epitope.

Reagents Provided

Bond Epitope Retrieval Solution 1 contains a citrate based buffer and surfactant. Total volume = 1 L.

pH 5.9–6.1 @ 25 °C

Dilution and Mixing

Bond Epitope Retrieval Solution 1 is ready to use. Do not dilute this reagent.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the bottle label. Return to 2–8 °C when not in use.

The signs indicating contamination and/or instability of Bond Epitope Retrieval Solution 1 are turbidity of the solution, odor development, and presence of precipitate.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.

Instructions for Use

Bond Epitope Retrieval Solution 1 should be poured into the bulk container marked "Epitope Retrieval Solution 1" located within the Bond Processing Module. This container can hold up to 1.3 L.

Product Specific Limitations

The appropriate time of heat retrieval may vary due to variation in tissue fixation and must be determined by the user.

Over retrieval of tissue sections may result in unwanted staining. Negative reagent controls should be used when optimizing retrieval conditions.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. To Vision

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