BOND FISH Kit

Catalog No: DS9636

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

EN FR IT DE ES PT
SV EL DA NL NO TR

Instructions for Use

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Catalog No: DS9636

Intended Use

This reagent is for *in vitro* diagnostic use.

The BOND FISH Kit enables the user to perform fluorescence *in situ* hybridization (FISH) on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system). It is intended for use with nucleic acid probes on formalin-fixed, paraffin-embedded (FFPE) tissue.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

The BOND FISH Kit consists of a formamide mixture recommended for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). This solution reduces non-specific hybridization of nucleic acid probes. The BOND FISH Kit enables FISH to be performed on the automated BOND system.

Reagents Provided

Post Hybridization Wash 2 Solution

Total volume = 18 mL, sufficient for 60 tests.

Dilution and Mixing

BOND FISH Kit is ready to use. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and *in situ* hybridization using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the container label.

There are no obvious signs that could indicate contamination and/or instability.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.

POST HYBRIDIZATION

H360D:

P201: Obtain special instructions before use.

WASH 2

May damage the unborn child.

P202: Do not handle until all safety precautions have been read and understood.

Contains Formamide (<50%).

P260: Do not breathe dust/fumes/gas/mist/vapours/spray.

GHS08: Health hazard.

P281: Use personal protective equipment as required.

Signal words: Danger.

P308+313: If exposed or concerned: Get medical advice/attention.

P314: Get medical advice/attention if you feel unwell.

Restricted to professional users.

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Measures to minimize microbial contamination of reagents need to be taken to prevent the occurrence of non-specific staining.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

BOND FISH Kit is developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) for use with BOND ancillary reagents and user selected FISH probes, including Kreatch XL Probes for BOND.

BOND FISH Kit should be utilized with BOND FISH Protocol D. This is designed as a post-hybridization wash protocol, to reduce non-specific hybridization of nucleic acid probes.

Please note that it is important for the user to set up appropriate pre-treatment and hybridization conditions. The test protocol may vary according to the probe selected by the user and by specimen type; please refer to the relevant Instructions for Use for further information.

Product Specific Limitations

Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. Protocol times may vary due to variation in tissue type, fixation and processing. In addition, BOND Enzyme concentration and incubation time may require optimization depending on probe, tissue type, processing and fixation conditions.

Performance Characteristics

The performance of BOND FISH Kit, DS9636, has been verified using a range of Kreatech XL Probes for BOND. The product is stable up to the expiry date printed on the product label.

Troubleshooting

Reference 3 may aid in remedial action.

Test samples should be complemented by the appropriate tissue controls.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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