# Novocastra™ Lyophilized Mouse Monoclonal Antibody CD57

**Product Code: NCL-NK1**

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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## Instructions for Use

Please read before using this product.

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## Check the integrity of the packaging before use.

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Kontroller, at pakken er ubeskadiget før brug.

# Novocastra™ Lyophilized Mouse Monoclonal Antibody CD57

## Product Code: NCL-NK1

### Intended Use

*For in vitro diagnostic use.*

NCL-NK1 is intended for the qualitative identification by light microscopy of CD57 molecules in paraffin sections. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Principle of Procedure

Immunohistochemical (IHC) staining techniques allow for the visualization of antigens via the sequential application of a specific antibody to the antigen (primary antibody), a secondary antibody to the primary antibody and an enzyme complex with a chromogenic substrate with interposed washing steps. The enzymatic activation of the chromogen results in a visible reaction product at the antigen site. The specimen may then be counterstained and coverslipped. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

### Clone

NK-1

### Immunogen

Human peripheral blood mononuclear cells.

### Specificity

Human CD57 molecule.

### Reagent Composition

NCL-NK1 is a lyophilized tissue culture supernatant containing 15 mM sodium azide as a preservative. The user is required to reconstitute the contents of the vial with the correct volume of sterile distilled water as indicated on the vial label.

### Ig class

IgM

### Total Protein Concentration Total Protein

1.0–8.0 g/L. Refer to vial label for lot specific total protein concentration.

### Antibody Concentration

Greater than or equal to 6.3 mg/L as determined by ELISA. Refer to vial label for batch specific Ig concentration.

### Recommendations On Use

Immunohistochemistry (see **D. Methodology**) on paraffin sections. Suggested dilution: 1:50 for 60 minutes at 25 °C. No pretreatment is recommended. This is provided as a guide and users should determine their own optimal working dilutions.

### Storage and Stability

Store unopened antibody at 2–8 °C. Under these conditions, there is no significant loss in product performance up to the expiry date indicated on the vial label. Do not use after expiration date indicated on the vial label. The reconstituted antibody is stable for at least two months when stored at 2–8 °C. For long term storage, it is recommended that aliquots of the reconstituted antibody are stored frozen at -20 °C (frost-free freezers are not recommended). Repeated freezing and thawing must be avoided. Prepare working dilutions on the day of use. Return to 2–8 °C immediately after use. Storage conditions other than those specified above must be verified by the user.

### Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

### Warnings and Precautions

This reagent has been prepared from the supernatant of cell culture. As it is a biological product, reasonable care should be taken when handling it.

The molarity of sodium azide in this reagent is 15 mM. A Material Safety Data Sheet (MSDS) is available upon request for sodium azide. Consult federal, state or local regulations for disposal of any potentially toxic components.

Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.<sup>1</sup> Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur.

Incubation times or temperatures, other than those specified, may give erroneous results. Any such changes must be validated by the user.

## Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures. Controls should be fresh autopsies/biopsy/surgical specimens, formalin-fixed, processed and paraffin wax-embedded as soon as possible in the same manner as the patient sample(s).

### Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques.

One positive tissue control should be included for each set of test conditions in each staining run.

A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.<sup>2</sup>

Recommended positive control tissue is tonsil.

If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

### Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target antigen by the primary antibody. Recommended negative control tissue is uterus.

Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user.

Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.<sup>3</sup> False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin (eg. liver, breast, brain, kidney) depending on the type of immunostain used. To differentiate endogenous enzyme activity or non-specific binding of enzymes from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate chromogen or enzyme complexes (avidin-biotin, streptavidin, labeled polymer) and substrate-chromogen, respectively. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

### Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

### Patient Tissue

Examine patient specimens stained with NCL-NK1 last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

## Results Expected

### Normal Tissues

Clone NK-1 detects the CD57 antigen on the membrane of lymphocytes in a variety of lymphoid and non-lymphoid tissues. Large granular lymphocytes in germinal centers of tonsil (16/16), spleen (9/9) and lymph node (1/1), were immunostained. Clone NK-1 also stained a number of different cell types in non-lymphoid tissues. In the tissues tested (n=99), staining was also observed in glandular epithelium of prostate, peripheral nerve and central nervous system, glandular epithelium of ovary, bronchiolar epithelium and in syncytiotrophoblasts. Variable staining was seen in the colloid of thyroid and some focal staining in kidney epithelium. Additionally, some staining was observed in capillaries of various tissues, subsets of pancreatic islet cells, hepatocytes and ovarian stromal cells.

### Abnormal Tissues

Clone NK-1 stained large scattered immunoreactive cells in T-cell lymphoma (5/5) and in B-cell lymphoma (5/5) whilst the neoplastic cells were negative. Additional staining was observed in a range of common tumors including colon adenocarcinoma, prostatic adenocarcinoma and ovarian carcinoma (9/26).

**NCL-NK1 is recommended for the identification of cells expressing CD57 in lymphoid and non-lymphoid tissues.**

### General Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artefacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.<sup>4</sup> Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Antibodies from Leica Biosystems Newcastle Ltd are for use, as indicated, on either frozen or paraffin-embedded sections with specific fixation requirements. Unexpected antigen expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

### Bibliography - General

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.

3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

**Amendments to Previous Issue**

Not applicable.

**Date of Issue**

23 June 2008 (NCL-NK1/CE/UK).

# Immunohistochemistry methodology for using Novocastra™ antibodies on paraffin-embedded tissue with ABC technique.

## Reagents required but not supplied

1. Standard solvents used in immunohistochemistry.
2. 0.5% v/v hydrogen peroxide.
3. 50 mM Tris-buffered saline (TBS) pH 7.6.
4. Antibody diluent - optimally diluted normal serum.
5. Normal sera from the species in which the secondary antibody is raised.
6. Secondary biotinylated antibody - prepare as recommended by manufacturer.
7. Avidin/Biotin Complex-Horseradish peroxidase (ABC-HRP) - prepare as recommended by manufacturer.
8. 3,3' Diaminobenzidine tetrahydrochloride (DAB) - prepare as recommended by manufacturer.
9. Hematoxylin counterstain - prepare as recommended by manufacturer.
10. Mounting medium - use as recommended by manufacturer.

## Equipment required but not supplied

1. Incubator set to 25 °C.
2. General immunohistochemistry laboratory equipment.

## Antigen retrieval solutions

Not applicable when no pretreatment is recommended.

## Methodology

**Prior to undertaking this methodology, users must be trained in immunohistochemical techniques.**

Customers should determine optimal dilutions for antibodies. Unless indicated, all steps are performed at room temperature (25 °C).

1. Cut and mount sections on slides coated with a suitable tissue adhesive.
2. De-paraffinize sections in xylene or xylene substitutes.
3. Re-hydrate through graded alcohols.
4. Neutralize endogenous peroxidase using 0.5% v/v hydrogen peroxide/methanol for 10 minutes.
5. Wash slides in running tap water.
6. Wash sections in TBS for 1 x 5 minutes with gentle rocking.
7. Cover sections with diluted normal serum for 10 minutes.
8. Incubate sections with optimally diluted primary antibody (see Recommendations on Use).
9. Wash in TBS buffer for 2 x 5 minutes with gentle rocking.
10. Incubate sections in appropriate biotinylated secondary antibody.
11. Wash in TBS buffer for 2 x 5 minutes with gentle rocking.
12. Incubate slides in ABC-HRP.
13. Wash in TBS buffer for 2 x 5 minutes with gentle rocking.
14. Incubate slides in DAB.
15. Rinse slides in water.
16. Counterstain with hematoxylin.
17. Dehydrate, clear and mount sections.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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## **Literatur - Allgemein**

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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## **Bibliografi - Allmän**

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
5. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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## Bibliografi - Generelt

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.

2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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