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# BOND™ Ready-To-Use Primary Antibody CD21 (2G9)

## Catalog No: PA0171

### Intended Use

*This reagent is for in vitro diagnostic use.*

CD21 (2G9) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human CD21 molecule in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). CD21 (2G9) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human CD21 molecule is achieved by first, allowing the binding of CD21 (2G9) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

### Reagents Provided

CD21 (2G9) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin® 950 as a preservative.

Total volume = 7 mL.

### Clone

2G9.

### Immunogen

Prokaryotic recombinant protein corresponding to the external domain of the human CD21 molecule.

### Specificity

Human CD21 molecule.

### Subclass

IgG2a.

### Total Protein Concentration

Approx 10 mg/mL.

### Antibody Concentration

Greater than or equal to 6.7 mg/L as determined by ELISA.

### Dilution and Mixing

CD21 (2G9) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of CD21 (2G9) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

### Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin® 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com).
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

CD21 (2G9) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The recommended staining protocol for CD21 (2G9) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 20 minutes.

## Results Expected

### Normal Tissues

Clone 2G9 stained the CD21 protein in follicular dendritic cells (FDCs) of lymphoid tissue and mature B cells. It strongly labelled FDCs in 29/29 cases of tonsil, with a dense meshed FDC network in the light zone and a more loosely arranged, weaker staining pattern of staining in the dark zone. The staining of mantle zone B cells was faint and much less consistent. The staining of the FDC network was also seen in 5/5 cases of spleen. FDC staining was also evident in 2/4 cases of thyroid, 2/2 cases of cervix, 3/4 cases of stomach and 3/4 cases of colon. The occasional mononuclear cell was labelled in other normal tissues tested. There was some weak staining of Leydig cells in testis no additional staining was observed in cases of uterus, placenta, prostate, ovary, pancreas, lung, liver, skin, breast, heart, small intestine, brain, pituitary and adrenal (n = 200).

### Tumor Tissues

With clone 2G9, in follicular lymphoma (n = 3) the atretic follicles exhibited strong, well-defined FDC networks. In the larger abnormal follicles a more diffuse and less intense staining of FDCs was observed. No expression of CD21 was observed in the neoplastic cells of T-cell lymphoma (n=6) other than a few residual follicles showing strong FDC staining. No staining of neoplastic cells was observed in non-Hodgkin's B-cell lymphoma (n = 14), lymphoblastic lymphoma (n = 2), Hodgkin's disease (n = 4) or a range of additional tumors (n=26).

**CD21 (2G9) is recommended for the immunodetection of normal and abnormal cell expression of CD21 in follicular dendritic cells.**

## Product Specific Limitations

CD21 (2G9) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kojima M, Nakamura S, Shimizu K, et al. Classical Hodgkin lymphoma occurring in clusters of nodal marginal zone B-lymphocytes in association with progressive transformation of germinal center. Pathology Research and Practice. 2003; 199:547–550.
5. Kojima M, Nakamura S, Motoori T, et al. Primary marginal zone B-cell lymphoma of the lymph node resembling plasmacytoma arising from a plasma cell variant of Castleman's disease. A clinicopathological and immunohistochemical study of seven patients. Acta Pathologica, Microbiologica et Immunologica Scandinavica (APMIS). 2002; 110:875–880.
6. Dupin N, Fisher C, Kellam P et al. Distribution of human herpesvirus-8 latently infected cells in Kaposi's sarcoma, multicentric Castleman's disease, and primary effusion lymphoma. Proceedings of the National Academy of Sciences USA. 1999; 96(8): 4546–4551.

ProClin® 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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## Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kojima M, Nakamura S, Shimizu K, et al. Classical Hodgkin lymphoma occurring in clusters of nodal marginal zone B-lymphocytes in association with progressive transformation of germinal center. Pathology Research and Practice. 2003; 199:547–550.
5. Kojima M, Nakamura S, Motoori T, et al. Primary marginal zone B-cell lymphoma of the lymph node resembling plasmacytoma arising from a plasma cell variant of Castleman's disease. A clinicopathological and immunohistochemical study of seven patients. Acta Pathologica, Microbiologica et Immunologica Scandinavica (APMIS). 2002; 110:875–880.
6. Dupin N, Fisher C, Kellam P et al. Distribution of human herpesvirus-8 latently infected cells in Kaposi's sarcoma, multicentric Castleman's disease, and primary effusion lymphoma. Proceedings of the National Academy of Sciences USA. 1999; 96(8):4546–4551.

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## Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kojima M, Nakamura S, Shimizu K, et al. Classical Hodgkin lymphoma occurring in clusters of nodal marginal zone B-lymphocytes in association with progressive transformation of germinal center. Pathology Research and Practice. 2003; 199:547-550.
5. Kojima M, Nakamura S, Motoori T, et al. Primary marginal zone B-cell lymphoma of the lymph node resembling plasmacytoma arising from a plasma cell variant of Castleman's disease. A clinicopathological and immunohistochemical study of seven patients. Acta Pathologica, Microbiologica et Immunologica Scandinavica (APMIS). 2002; 110:875-880.
6. Dupin N, Fisher C, Kellam P et al. Distribution of human herpesvirus-8 latently infected cells in Kaposi's sarcoma, multicentric Castleman's disease, and primary effusion lymphoma. Proceedings of the National Academy of Sciences USA. 1999; 96(8):4546-4551.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kojima M, Nakamura S, Shimizu K, et al. Classical Hodgkin lymphoma occurring in clusters of nodal marginal zone B-lymphocytes in association with progressive transformation of germinal center. Pathology Research and Practice. 2003; 199:547-550.
5. Kojima M, Nakamura S, Motoori T, et al. Primary marginal zone B-cell lymphoma of the lymph node resembling plasmacytoma arising from a plasma cell variant of Castleman's disease. A clinicopathological and immunohistochemical study of seven patients. Acta Pathologica, Microbiologica et Immunologica Scandinavica (APMIS). 2002; 110:875-880.
6. Dupin N, Fisher C, Kellam P et al. Distribution of human herpesvirus-8 latently infected cells in Kaposi's sarcoma, multicentric Castleman's disease, and primary effusion lymphoma. Proceedings of the National Academy of Sciences USA. 1999; 96(8): 4546-4551.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kojima M, Nakamura S, Shimizu K, et al. Classical Hodgkin lymphoma occurring in clusters of nodal marginal zone B-lymphocytes in association with progressive transformation of germinal center. Pathology Research and Practice. 2003; 199:547–550.

5. Kojima M, Nakamura S, Motoori T, et al. Primary marginal zone B-cell lymphoma of the lymph node resembling plasmacytoma arising from a plasma cell variant of Castleman's disease. A clinicopathological and immunohistochemical study of seven patients. *Acta Pathologica, Microbiologica et Immunologica Scandinavica (APMIS)*. 2002; 110:875–880.
6. Dupin N, Fisher C, Kellam P et al. Distribution of human herpesvirus-8 latently infected cells in Kaposi's sarcoma, multicentric Castleman's disease, and primary effusion lymphoma. *Proceedings of the National Academy of Sciences USA*. 1999; 96(8): 4546–4551.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd