Bond™ Ready-to-Use Primary Antibody Lambda Light Chain (SHL53)

Catalog No: PA0570

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

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Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Kontroller, at pakken er ubeskadiget før brug.

Bond™ Ready-To-Use Primary Antibody Lambda Light Chain (SHL53)

Catalog No: PA0570

Intended Use

This reagent is for in vitro diagnostic use.

Lambda Light Chain (SHL53) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human lambda light chains in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated Bond system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation). The demonstration of human lambda light chains is achieved by first allowing the binding of Lambda Light Chain (SHL53) to the section and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated Bond system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Lambda Light Chain (SHL53) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

SHL53.

Immunogen

Prokaryotic recombinant protein corresponding to 105 amino acids of the human lambda light chain molecule.

Specificity

Human lambda light chain.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 30 mg/L as determined by ELISA.

Dilution and Mixing

Lambda Light Chain (SHL53) primary antibody is optimally diluted for use on the automated Bond system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the automated Bond system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Lambda Light Chain (SHL53) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Microsystems, or alternatively, visit the Leica Microsystems Web site, www.leica-microsystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.

- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Lambda Light Chain (SHL53) primary antibody was developed for use on the automated Bond system in combination with Bond Polymer Refine Detection. The recommended staining protocol for Lambda Light Chain (SHL53) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using Bond Epitope Retrieval Solution 1 for 30 minutes.

Results Expected

Normal Tissues

Clone SHL53 detected lambda light chains in serum and a variety of tissues/cells types. In lymphoid tissues, lambda positive plasma cells stained strongly, with weaker staining of lambda positive B-lymphocytes. (Total number of cases stained = 95).

Tumor Tissues

Clone SHL53 detected the lambda light chain protein in 57/72 hematological malignancies evaluated, including 14/14 diffuse large B cell lymphomas, 10/13 Hodgkin's lymphomas, 7/13 chronic lymphocytic leukemias, 6/6 mantle cell lymphomas, 5/6 follicular lymphomas, 4/6 peripheral T cell lymphomas, 2/2 anaplastic large cell lymphomas, 2/2 MALTomas, 4/4 unspecified lymphomas, 1/1 unspecified lymphoma of the spleen, 1/1 unspecified malignant lymphomas, 1/1 Burkitt's lymphomas, 1/1 acute lymphocytic leukemias, 1/1 marginal zone lymphomas, 0/1 thymomas, 0/1 NK/T cell lymphomas and 0/1 angioimmunoblastic lymphomas. Except for infiltrating lymphocytes, the lambda light chain was also detected in 30/40 non-lymphoid malignancies including 10/10 squamous cell carcinomas, 4/4 liver tumors, 2/4 ovarian tumors, 2/3 colorectal tumors, 2/2 breast tumors, 2/2 stomach tumors, 2/2 pancreatic tumors, 2/2 kidney tumors, 2/2 prostate tumors, 1/2 lung tumors and 1/2 lymph node metastatic carcinomas. No staining was detected in thyroid tumors (0/3) and brain tumors (0/2). (Total number of cases =112).

Lambda Light Chain (SHL53) is recommended for the detection of lambda light chains in normal and neoplastic tissues.

Product Specific Limitations

Lambda Light Chain (SHL53) has been optimized at Leica Microsystems for use with Bond Polymer Refine Detection and Bond ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Microsystems to report unusual staining.

Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 est une marque commerciale de Supelco, membre du groupe Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 è un marchio di fabbrica di Supelco, società del gruppo Sigma-Aldrich.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 ist eine Marke von Supelco, einem Unternehmen der Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 es una marca registrada de Supelco, parte de Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 é uma marca registada de Supelco, parte da Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 är ett varumärke tillhörande Supelco som ingår i Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. Το ProClim™ 950 είναι εμπορικό σήμα της Supelco, τμήμα της Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.

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17 September 2010