Bond™ Ready-to-Use Primary Antibody Human Follicle Stimulating Hormone (INN-hFSH-60)

Catalog No: PA0693

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Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Catalog No: PA0693

Intended Use

This reagent is for in vitro diagnostic use.

Human Follicle Stimulating Hormone (INN-hFSH-60) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human follicle stimulating hormone in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated Bond™ system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation). Human Follicle Stimulating Hormone (INN-hFSH-60) primary antibody is a ready to use product that has been specifically optimized for use with Bond Polymer Refine Detection. The demonstration of human follicle stimulating hormone is achieved by first, allowing the binding of Human Follicle Stimulating Hormone (INN-hFSH-60) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated Bond system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Human Follicle Stimulating Hormone (INN-hFSH-60) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

INN-hFSH-60.

Immunogen

Human follicle stimulating hormone.

Specificity

The beta 2 epitope of the human follicle stimulating hormone, beta subunit. Not crossreactive with human thyroid stimulating hormone, human luteinizing hormone, human chorionic gonadotrophin, alpha human chorionic gonadotrophin and alpha human follicle stimulating hormone.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.08 mg/L.

Dilution and Mixing

Human Follicle Stimulating Hormone (INN-hFSH-60) primary antibody is optimally diluted for use on the Bond system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Human Follicle Stimulating Hormone (INN-hFSH-60) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Human Follicle Stimulating Hormone (INN-hFSH-60) primary antibody was developed for use on the automated Bond system in combination with Bond Polymer Refine Detection. The recommended staining protocol for Human Follicle Stimulating Hormone (INN-hFSH-60) primary antibody is IHC Protocol F. Enzyme pretreatment is recommended using Bond Enzyme 1 for 10 minutes.

Results Expected

Normal Tissues

Clone INN-hFSH-60 detected human follicle stimulating hormone in gonadotropic cells of the anterior pituitary gland (Total number of normal tissues = 66).

Tumor Tissues

Clone INN-hFSH-60 gave no tumor specific staining in the cases evaluated, including squamous cell carcinomas, adenocarcinomas, sarcomas, melanomas and lymphomas (Total number of cases stained = 44).

Human Follicle Stimulating Hormone (INN-hFSH-60) is recommended for the detection of human follicle stimulating hormone by qualitative immunohistochemistry on normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

Product Specific Limitations

Human Follicle Stimulating Hormone (INN-hFSH-60) has been optimized at Leica Biosystems for use with Bond Polymer Refine Detection and Bond ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Holm R, Nesland J M, Attramadal A, et al. Expression of Leptin Receptor Isoforms and Effects of Leptin on the Proliferation and Hormonal Secretion in Human Pituitary Adenomas. Journal of Pathology.1989; 158:213–217.
5. Jautzke G. Simultaneous production of the alpha-subunit of glycoprotein hormones and other hormones in pituitary adenomas. Pathology – Research and Practice. 1988; 183(5):601–605.
6. Hitchcock E and Morris C S. Thyroid stimulating and gonadotrophic hormones in pituitary adenomas without clinical or serological abnormality. Histochemical Journal. 1986; 18(6):317–320.

ProClin® 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Holm R, Nesland J M, Attramadala A, et al. Expression of Leptin Receptor Isoforms and Effects of Leptin on the Proliferation and Hormonal Secretion in Human Pituitary Adenomas. Journal of Pathology. 1989; 158:213–217.
5. Jautzke G. Simultaneous production of the alpha-subunit of glycoprotein hormones and other hormones in pituitary adenomas. Pathology – Research and Practice. 1988; 183(5):601–605.
6. Hitchcock E and Morris C S. Thyroid stimulating and gonadotrophic hormones in pituitary adenomas without clinical or serological abnormality. Histochemical Journal. 1986; 18(6):317–320.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Holm R, Nesland J M, Attramadal A, et al. Expression of Leptin Receptor Isoforms and Effects of Leptin on the Proliferation and Hormonal Secretion in Human Pituitary Adenomas. Journal of Pathology. 1989; 158:213–217.
5. Jautzke G. Simultaneous production of the alpha-subunit of glycoprotein hormones and other hormones in pituitary adenomas. Pathology – Research and Practice. 1988; 183(5):601–605.
6. Hitchcock E and Morris C S. Thyroid stimulating and gonadotrophic hormones in pituitary adenomas without clinical or serological abnormality. Histochemical Journal. 1986; 18(6):317–320.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCLCS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Holm R, Nesland J M, Attramadala A, et al. Expression of Leptin Receptor Isoforms and Effects of Leptin on the Proliferation and Hormonal Secretion in Human Pituitary Adenomas. Journal of Pathology. 1989; 158:213–217.
5. Jautzke G. Simultaneous production of the alpha-subunit of glycoprotein hormones and other hormones in pituitary adenomas. Pathology – Research and Practice. 1988; 183(5):601–605.
6. Hitchcock E and Morris C S. Thyroid stimulating and gonadotrophic hormones in pituitary adenomas without clinical or serological abnormality. Histochemical Journal. 1986; 18(6):317–320.

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ProClin[®] 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.