# Bond™ Ready-to-Use Primary Antibody Cytokeratin 20 (PW31)

**Catalog No: PA0918**

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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## **Instructions for Use**

Please read before using this product.

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## **Check the integrity of the packaging before use.**

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Kontroller, at pakken er ubeskadiget før brug.

# Bond™ Ready-To-Use Primary Antibody Cytokeratin 20 (PW31)

## Catalog No: PA0918

### Intended Use

*This reagent is for in vitro diagnostic use.*

Cytokeratin 20 (PW31) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human cytokeratin 20 intermediate filament protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated Bond™ system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation). Cytokeratin 20 (PW31) primary antibody is a ready to use product that has been specifically optimized for use with Bond Polymer Refine Detection. The demonstration of human cytokeratin 20 intermediate filament protein is achieved by first, allowing the binding of Cytokeratin 20 (PW31) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated Bond system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

### Reagents Provided

Cytokeratin 20 (PW31) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

### Clone

PW31.

### Immunogen

Prokaryotic recombinant protein corresponding to a 70 amino acid component of the N-terminal region of the cytokeratin 20 intermediate filament protein.

### Specificity

Human cytokeratin 20 intermediate filament protein.

### Subclass

IgG1.

### Total Protein Concentration

Approx 10 mg/mL.

### Antibody Concentration

Greater than or equal to 1.5 mg/L as determined by ELISA.

### Dilution and Mixing

Cytokeratin 20 (PW31) primary antibody is optimally diluted for use on the Bond system. Reconstitution, mixing, dilution or titration of this reagent is not required.

### Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond system.

### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Cytokeratin 20 (PW31) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

## Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com).
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

Cytokeratin 20 (PW31) primary antibody was developed for use on the automated Bond system in combination with Bond Polymer Refine Detection. The recommended staining protocol for Cytokeratin 20 (PW31) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using Bond Epitope Retrieval Solution 1 for 20 minutes.

## Results Expected

### Normal Tissues

Clone PW31 detects cytokeratin 20 in the cytoplasm of normal gastric, small and large bowel epithelium, urothelium.

### Tumor Tissues

Clone PW31 stained carcinomas of the colon (60/65) and bladder (1/2) and Merkel cell tumours. Focal immunostaining was seen in lung squamous cell carcinomas (3/7). It did not stain a range of other tumours (0/99), including breast carcinomas, renal clear cell carcinoma, prostate adenocarcinomas, pancreatic adenocarcinoma, endometrial carcinoma, thyroid carcinomas, lung small cell carcinomas and adenocarcinomas, melanomas and sarcomas.

Cytokeratin 20 (PW31) is recommended for use as part of an antibody panel in the characterization of tumors originating from urothelium, intestinal epithelium and Merkel cells.

## Product Specific Limitations

Cytokeratin 20 (PW31) has been optimized at Leica Biosystems for use with Bond Polymer Refine Detection and Bond ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leech SN, Kolar AJO, Barrett PD et al. Merkel cell carcinoma can be distinguished from metastatic small cell carcinoma using antibodies to cytokeratin 20 and thyroid transcription factor 1. Journal of Clinical Pathology. 2001; 54:727–729.
5. Campbell F and Herrington CS. Application of cytokeratin 7 and 20 immunohistochemistry to diagnostic pathology. Current Diagnostic Pathology. 2001; 7:113–122.

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## Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Référence à commander : M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4ème édition. Churchill Livingstone, New York. 1996.
4. Leech SN, Kolar AJO, Barrett PD et al. Merkel cell carcinoma can be distinguished from metastatic small cell carcinoma using antibodies to cytokeratin 20 and thyroid transcription factor 1. Journal of Clinical Pathology 2001;54:727-729.
5. Campbell F and Herrington CS. Application of cytokeratin 7 and 20 immunohistochemistry to diagnostic pathology. Current Diagnostic Pathology 2001;7:113-122.

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## Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 28. Februar 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD und Stevens A. Theory and Practice of Histological Techniques. 4. Auflage. Churchill Livingstone, New York. 1996.
4. Leech SN, Kolar AJO, Barrett PD et al. Merkel cell carcinoma can be distinguished from metastatic small cell carcinoma using antibodies to cytokeratin 20 and thyroid transcription factor 1. Journal of Clinical Pathology. 2001; 54:727–729.
5. Campbell F und Herrington CS. Application of cytokeratin 7 and 20 immunohistochemistry to diagnostic pathology. Current Diagnostic Pathology. 2001; 7:113–122.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leech SN, Kolar AJO, Barrett PD et al. Merkel cell carcinoma can be distinguished from metastatic small cell carcinoma using antibodies to cytokeratin 20 and thyroid transcription factor 1. Journal of Clinical Pathology. 2001; 54:727–729.
5. Campbell F and Herrington CS. Application of cytokeratin 7 and 20 immunohistochemistry to diagnostic pathology. Current Diagnostic Pathology. 2001; 7:113–122.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leech SN, Kolar AJO, Barrett PD et al. Merkel cell carcinoma can be distinguished from metastatic small cell carcinoma using antibodies to cytokeratin 20 and thyroid transcription factor 1. Journal of Clinical Pathology. 2001; 54:727-729.
5. Campbell F and Herrington CS. Application of cytokeratin 7 and 20 immunohistochemistry to diagnostic pathology. Current Diagnostic Pathology. 2001; 7:113-122.

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## Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leech SN, Kolar AJO, Barrett PD et al. Merkel cell carcinoma can be distinguished from metastatic small cell carcinoma using antibodies to cytotkeratin 20 and thyroid transcription factor 1. Journal of Clinical Pathology. 2001; 54:727–729.
5. Campbell F and Herrington CS. Application of cytotkeratin 7 and 20 immunohistochemistry to diagnostic pathology. Current Diagnostic Pathology. 2001; 7:113–122.

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27 juni 2014