Bond™ Ready-to-Use Primary Antibody CD43 (MT1)

Catalog No: PA0938

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Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Catalog No: PA0938

Intended Use

This reagent is for in vitro diagnostic use.

CD43 (MT1) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human CD43 antigen in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated Bond-max™ system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation). CD43 (MT1) primary antibody is a ready to use product that has been specifically optimized for use with Bond Polymer Refine Detection. The demonstration of human CD43 antigen is achieved by first, allowing the binding of CD43 (MT1) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated Bond-max system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

CD43 (MT1) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

MT1.

Immunogen

Human lymphocytes.

Specificity

Human CD43 antigen.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.3 mg/L as determined by ELISA.

Dilution and Mixing

CD43 (MT1) primary antibody is optimally diluted for use on the Bond-max system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond-max system.

Storage and Stability

Store at 2-8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of CD43 (MT1) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2-8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

CD43 (MT1) primary antibody was developed for use on the automated Bond-max system in combination with Bond Polymer Refine Detection. The recommended staining protocol for CD43 (MT1) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using Bond Epitope Retrieval Solution 2 for 20 minutes.

Results Expected

Normal Tissues

Clone MT1 detected the CD43 antigen on the membrane and in the cytoplasm of T cells and cells of myeloid origin. Staining was not seen in non-hemopoietic tissue elements. (n=85).

Tumor Tissues

Clone MT1 stained 21/43 lymphomas, including 5/14 diffuse large B cell lymphomas, 1/6 follicle center lymphomas, 0/5 Hodgkin's disease, 2/2 mantle cell lymphomas, 1/3 MALTomas, 6/6 peripheral T cell lymphomas, 1/2 anaplastic large cell lymphomas, 1/1 Burkitt's lymphoma, 1/1 T lymphoblastic lymphoma, 1/1 angioimmunoblastic T cell lymphoma, 1/1 NK/T cell lymphoma, and 1/1 T cell rich B cell lymphoma. Non-lymphoid tumors (n=41) were negative, except for infiltrating T cells.

CD43 (MT1) is recommended for use as part of a panel of antibodies for immunophenotyping T and B cell lymphomas and the identification of myeloid disorders.

Product Specific Limitations

CD43 (MT1) has been optimized at Leica Biosystems for use with Bond Polymer Refine Detection and Bond ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Attygalle AD, Liu H, Shirali S, et al. Atypical marginal zone hyperplasia of mucosa-associated lymphoid tissue: a reactive condition of childhood showing immunoglobulin lambda light-chain restriction. Blood. 2004; 104(10):3343-3348.
5. Itoh T, Shimizu M, Kitami K, et al. Primary extranodal marginal zone B-cell lymphoma of the mucosa-associated lymphoid tissue type in the CNS. Neuropathology. 2001; 21:174-180.
6. Bulmer JN, Jones RK and Searle RF. Intraepithelial leukocytes in endometriosis and adenomyosis: comparison of eutopic and ectopic endometrium with normal endometrium. Human Reproduction. 1998; 13(10):2910-2915.

ProClin™ 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD och Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Attygalle AD, Liu H, Shirali S, et al. Atypical marginal zone hyperplasia of mucosa-associated lymphoid tissue: a reactive condition of childhood showing immunoglobulin lambda light-chain restriction. Blood. 2004; 104(10):3343-3348.
5. Itoh T, Shimizu M, Kitami K, et al. Primary extranodal marginal zone B-cell lymphoma of the mucosa-associated lymphoid tissue type in the CNS. Neuropathology. 2001; 21:174-180.
6. Bulmer JN, Jones RK och Searle RF. Intraepithelial leukocytes in endometriosis and adenomyosis: comparison of eutopic and ectopic endometrium with normal endometrium. Human Reproduction. 1998; 13(10):2910-2915.

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Leica Biosystems og logodesign er varemærker tilhørende Leica Biosystems A.G. og anvendes under licens. © 2007 Leica Biosystems Newcastle Ltd. Alle rettigheder forbeholdes.
ProClin™ 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.
Leica Biosystems Newcastle Ltd