

# Novocastra™ Lyophilized Mouse Monoclonal Antibody Plasminogen Activator Inhibitor (Type 1)

## Product Code: NCL-PAI-1

<b>Intended Use</b>	<b>For In Vitro Diagnostic Use:</b> This product is intended for qualitative immunohistochemistry with normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.
<b>Specificity</b>	Human plasminogen activator inhibitor type 1 (PAI-1).
<b>Clone</b>	TJA6
<b>Ig Class</b>	IgG2b, kappa
<b>Antigen Used for Immunizations</b>	Prokaryotic recombinant protein corresponding to a 250 amino acid portion of the N-terminus of the human PAI-1 molecule.
<b>Hybridoma Partner</b>	Mouse myeloma (p3-NS1-Ag4-1).
<b>Preparation</b>	Lyophilized tissue culture supernatant containing 15 mM sodium azide. Reconstitute with 1 mL or 0.1 mL of sterile distilled water as indicated on vial label.
<b>Effective on Frozen Tissue</b>	Not fully evaluated.
<b>Effective on Paraffin Wax Embedded Tissue</b>	Yes
<b>Recommendations on Use</b>	Immunohistochemistry: Typical working dilution 1:20–1:40. 60 minutes primary antibody incubation at 25 °C. Standard ABC technique. Western Blotting: Not recommended.
<b>Positive Controls</b>	Immunohistochemistry: Tonsil.
<b>Staining Pattern</b>	Cytoplasmic.
<b>Storage and Stability</b>	Store unopened lyophilized antibody at 4 °C. Under these conditions, there is no significant loss in product performance up to the expiry date indicated on the vial label. The reconstituted antibody is stable for at least two months when stored at 4 °C. For long term storage, it is recommended that aliquots of the antibody are frozen at -20 °C (frost-free freezers are not recommended). Repeated freezing and thawing must be avoided. Prepare working dilutions on the day of use.
<b>General Overview</b>	Plasminogen activator inhibitor type 1 (PAI-1) is a 48 kD protein which inhibits the conversion of plasminogen to plasmin. It is the principal inhibitor of the plasminogen activators t-PA and u-PA. PAI-1 is structurally related to the serine protease inhibitor (serpin) superfamily. The serpins are known to undergo a conformational rearrangement upon cleavage of the reactive central peptide bond (P1-P1) and it is this conformational difference between the active and cleaved forms which determines their reactivity.
<b>General References</b>	Dublin E, Hanby A, Patel N K, et al. <i>American Journal of Pathology</i> . 157 (4): 1219–1227 (2000). de Witte J H, Sweep C G, Klijn J G, et al. <i>British Journal of Cancer</i> . 80 (1–2): 286–294 (1999). Hartbeck N, Thomssen C, Berger U, et al. <i>Breast Cancer Research Treatment</i> . 54 (2): 147–157 (1999). Arai Y, Kubota T, Nakagawa T, et al. <i>Acta. Neurochir. (Wien)</i> . 140 (4): 377–385 (1998). Strojan P, Budihna M, Smid L, et al. <i>European Journal of Cancer</i> . 34 (8): 1193–1197 (1998). Sweep C G, Geurts-Moespot J, Grebenschikov N, et al. <i>British Journal of Cancer</i> . 78 (11): 1434–1441 (1998). Speiser P, Mayerhofer K, Kucera E, et al. <i>Anticancer Research</i> . 17 (1B): 679–683 (1997). Ito H, Yonemura Y, Fujita H, et al. <i>Virchows Arch</i> . 427 (5): 487–496 (1996). Zellinger R, Eder S, Schneeberger C, et al. <i>Anticancer Research</i> . 16 (1): 449–453 (1996). Sancho E, Declerck P J, Price N C, et al. <i>Biochemistry</i> . 34 (3): 1064–1069 (1995). Keijer J, Linders M, van Zonneveld A-J, et al. <i>Blood</i> . 78 (2): 401–409 (1991). Munch M, Heegaard C, Jensen P H, et al. <i>FEBS</i> . 295 (1,2,3): 102–106 (1991).

