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Check the integrity of the packaging before use.

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Intended Use

For *in vitro* diagnostic use.

Novolink™ Polymer Detection Systems are used for the visualization of mouse IgG, mouse IgM and rabbit IgG primary antibodies. Novolink™ Polymer and Novolink™ DAB (Polymer) contain component reagents of these systems and are intended for use in the procedure described below. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Principle of Procedure

The first immunohistoperoxidase technique was reported by Nakane and Pierce.¹ Since then many developments have occurred, leading to increased sensitivity over earlier techniques. A recent development has been the use of polymeric labeling. This technology has been applied to both primary antibodies² and detection systems. The Novolink™ Polymer Detection Systems utilize a novel controlled polymerization technology to prepare polymeric HRP-linker antibody conjugates. Therefore, the problem of non-specific staining that can occur with Streptavidin/Biotin detection systems due to endogenous biotin does not occur.

These products are used in an immunohistochemical (IHC) procedure, which allows the qualitative identification by light microscopy of antigens in sections of formalin-fixed, paraffin-embedded tissue, via sequential steps with interposed washing steps. If required by the primary antibody, sections are subjected to epitope retrieval prior to staining. Endogenous peroxidase activity is neutralized using the Peroxidase Block. This is followed by application of the Novocastra™ Protein Block to reduce non-specific binding of primary and polymer. The section is subsequently incubated with optimally diluted primary antibody. Post Primary (Rabbit anti Mouse IgG) is then used to detect mouse antibodies. The Novolink™ Polymer recognizes rabbit immunoglobulins, it detects the post primary and any tissue-bound rabbit primary antibodies. Sections are further incubated with the substrate/chromogen, 3,3' - diaminobenzidine (DAB), prepared from DAB Chromogen and Novolink™ DAB Substrate Buffer (Polymer), as described below. Reaction with the peroxidase produces a visible brown precipitate at the antigen site. Sections are counterstained with Hematoxylin and coverslipped. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

Reagents Provided

Details of which reagents from the following list are supplied in each product are given in the table below.

1. Peroxidase Block. 3–4%(v/v) Hydrogen peroxide.
2. Protein Block. 0.4% Casein in phosphate-buffered saline, with stabilizers, surfactant, and 0.2% Bronidox L as a preservative.
3. Post Primary. Rabbit anti mouse IgG (<10 µg/mL) in 10% (v/v) animal serum in tris-buffered saline/0.09% ProClin™ 950.
4. Novolink™ Polymer. Anti-rabbit Poly-HRP-IgG (<25µg/mL) containing 10% (v/v) animal serum in tris-buffered saline/0.09% ProClin™ 950.
5. DAB Chromogen. 1.74% w/v 3,3' - diaminobenzidine, in a stabilizer solution.
6. Novolink™ DAB Substrate Buffer (Polymer). Buffered solution containing ≤0.1% hydrogen peroxide and preservative.
7. Hematoxylin. <0.1% Hematoxylin.

Reagents	Product Number	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagents	Product Number	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

Reconstitution, Mixing, Dilution, Titration

The Peroxidase Block, Protein Block, Post Primary, Novolink™ Polymer and Hematoxylin are prediluted. Reconstitution, mixing, dilution, or titration of these reagents is not recommended. Further dilution may result in loss of antigen staining. The user must validate any such change.

The DAB Chromogen requires dilution to 1/20 in Novolink™ DAB Substrate Buffer (Polymer) prior to use. Further dilution may result in loss of antigen staining. The user must validate any such change.

Storage and Stability

Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the product label. Storage conditions other than those specified must be verified by the user. There are no obvious signs to indicate instability of this product therefore positive and negative controls should be run simultaneously with patient samples.

Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

Warnings and Precautions

DAB Chromogen

Contains <10% Biphenyl-3,3', 4,4'-Tetrahydroamine.

GHS08: Health hazard.

Signal words: Danger.

H341: Suspected of causing genetic defects.

H350: May cause cancer.

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P281: Use personal protective equipment as required.

P308+313: IF exposed or concerned: Get medical advice/attention.

P501: Dispose of contents/container to hazardous or special waste collection point.

A Material Safety Data Sheet is available upon request or available from www.LeicaBiosystems.com

For professional users.

Do not mix reagents from different detection systems.

Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.⁴

Never pipette reagents by mouth and avoid contacting the skin and mucus membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.

Consult Federal, State or local regulations for disposal of any potentially toxic components.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur. Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Procedure

A. Reagents required but not supplied

1. Standard solvents used in immunohistochemistry.
2. 50mM tris-buffered saline (TBS) pH7.6.
3. Antigen retrieval solution(s).
4. Enzyme retrieval solution(s).
5. Antibody diluent.
6. Primary antibody.
7. Mounting medium.

B. Equipment required but not supplied

1. Equipment required for antigen retrieval, if recommended for the primary antibody.
2. General immunohistochemistry laboratory equipment.

C. Methodology

Prior to undertaking this methodology, users must be trained in immunohistochemical techniques.

All steps must be followed as directed or performance may be impaired.

The combination of the primary antibody, its dilution, together with the detection system should be validated by the user on a series of known positive and negative controls.

Unless indicated, all steps are performed at room temperature (25 °C).

For use on frozen tissue, cut sections and fix according to recommendations for primary antibody, commence at step 11.

1. Cut and mount sections on slides coated with a suitable tissue adhesive.
2. De-paraffinize sections in xylene or xylene substitutes.
3. Re-hydrate through graded alcohols.
4. Wash slides in running tap water.
5. Perform antigen retrieval as required (see **Recommendations for Use** for primary antibody).
6. Wash slides in de-ionized water.
7. Neutralize endogenous peroxidase using Peroxidase Block for 5 minutes.
8. Wash in TBS for 2 x 5 minutes.

9. Incubate with Protein Block for 5 minutes.
10. Wash in TBS for 2 x 5 minutes.
11. Incubate with optimally diluted primary antibody (see **Recommendations for Use** for primary antibody).
12. Wash in TBS for 2 x 5 minutes.
13. Incubate with Post Primary for 30 minutes.
14. Wash in TBS for 2 x 5 minutes.
15. Incubate with Novolink™ Polymer for 30 minutes.
16. Wash in TBS for 2 x 5 minutes with gentle rocking.
17. Develop peroxidase activity with DAB working solution (see DAB Working Solution) for 5 minutes.
18. Rinse slides in water.
19. Counterstain with Hematoxylin.
20. Rinse slides in water for 5 minutes.
21. Dehydrate, clear and mount sections.

DAB Working Solution

Add 50µl of DAB Chromogen to 1ml of Novolink™ DAB Substrate Buffer (Polymer). Use within six hours of preparation.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures. Controls should be fresh autopsy/biopsy/surgical specimens, formalin-fixed, processed and paraffin wax-embedded as soon as possible in the same manner as the patient sample(s).

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques. One positive tissue control should be included for each set of test conditions/primary antibody in each staining run. A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.⁵ For recommended positive control tissue see primary antibody Recommendations for Use. If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target antigen by the primary antibody. For recommended negative control tissue see primary antibody Instructions for Use. Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user. Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.⁶ False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin⁷ (eg. liver, breast, brain, kidney). To differentiate endogenous enzyme activity or non-specific binding from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate-chromogen, with labeled polymer and substrate-chromogen or with Post Primary, labeled polymer and substrate-chromogen. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

Patient Tissue

Examine patient specimens last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.⁸

Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Novolink™ Polymer Detection Systems and their components are for use on paraffin-embedded sections with specific fixation requirements. Unexpected antigen expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

Performance Characteristics

The performance of Novolink™ Polymer Detection Systems, Novolink™ Polymer and Novolink™ DAB (Polymer) have been validated using a range of Novocastra mouse IgG, mouse IgM* and rabbit IgG primary antibodies.

*Weak staining may be seen with some antibodies of IgM isotype.

These products are stable up to the expiry date(s) printed on the product label.

Bibliography

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of antigens. *Journal of Histochemistry and Cytochemistry*. 1967; 14:929–931.
2. Tsutsumi Y, Serizawa A and Kawai K. Enhanced polymer one-step staining (EPOS) for proliferating cell nuclear antigen and Ki-67 antigen-applications to intraoperative frozen diagnosis. *Pathology International*. 1995; 45(2):108–115.
3. Shan-Rong Shi, Guo J Cote R J et al . *Applied Immunohistochemistry & Molecular Morphology*. 1999; 7:201-208.
4. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
5. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
6. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
7. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
8. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Réactifs	Référence du produit	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Réactifs	Référence du produit	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Bibliographie

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of antigens. *Journal of Histochemistry and Cytochemistry*. 1967; 14:929–931.
2. Tsutsumi Y, Serizawa A and Kawai K. Enhanced polymer one-step staining (EPOS) for proliferating cell nuclear antigen and Ki-67 antigen-applications to intraoperative frozen diagnosis. *Pathology International*. 1995; 45(2):108–115.
3. Shan-Rong Shi, Guo J Cote R J et al . *Applied Immunohistochemistry & Molecular morphology*. 1999; 7:201-208.
4. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
5. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
6. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
7. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
8. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reagenti	Cod. prodotto	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagenti	Cod. prodotto	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Reagenzien	Produkt Nr.	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagenzien	Produkt Nr.	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Literatur

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of antigens. *Journal of Histochemistry and Cytochemistry*. 1967; 14:929–931.
2. Tsutsumi Y, Serizawa A and Kawai K. Enhanced polymer one-step staining (EPOS) for proliferating cell nuclear antigen and Ki-67 antigen-applications to intraoperative frozen diagnosis. *Pathology International*. 1995; 45(2):108–115.
3. Shan-Rong Shi, Guo J Cote R J et al . *Applied Immunohistochemistry & Molecular morphology*. 1999; 7:201-208.
4. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
5. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc. Philadelphia.
6. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
7. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
8. Omata M, Liew CT, Ashcavi M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reactivos	Referencia	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reactivos	Referencia	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of antigens. *Journal of Histochemistry and Cytochemistry*. 1967; 14:929–931.
2. Tsutsumi Y, Serizawa A and Kawai K. Enhanced polymer one-step staining (EPOS) for proliferating cell nuclear antigen and Ki-67 antigen-applications to intraoperative frozen diagnosis. *Pathology International*. 1995; 45(2):108–115.
3. Shan-Rong Shi, Guo J Cote R J et al . *Applied Immunohistochemistry & Molecular morphology*. 1999; 7:201-208.
4. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
5. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc. Philadelphia.
6. P6. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
7. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
8. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reagentes	N.º do produto	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagentes	N.º do produto	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Reagens	Produktnr.	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagens	Produktnr.	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Reagenser	Produkt Nr.	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagenser	Produkt Nr.	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Αντιδραστήρια	Κωδικός είδους:	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Αντιδραστήρια	Κωδικός είδους:	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Reagentia	Productnr.	Novolink™ Max Polymer Detection System (1250 tests) RE7280-K	Novolink™ Polymer Detection System (500 tests) RE7150-K	Novolink™ Polymer Detection System (250 tests) RE7140-K	Novolink™ Min Polymer Detection System (50 tests) RE7290-K
Peroxidase Block	RE7101		2 x 25ml	1 x 25ml	
Protein Block	RE7102		2 x 25ml	1 x 25ml	
Post Primary	RE7111		2 x 25ml	1 x 25ml	
Novolink™ Polymer	RE7112		2 x 25ml	1 x 25ml	
DAB Chromogen	RE7105		1 x 3ml	1 x 3ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7143		2 x 30ml	1 x 30ml	
Hematoxylin	RE7107		2 x 25ml	1 x 25ml	
Peroxidase Block	RE7157	1 x 125ml			
Protein Block	RE7158	1 x 125ml			
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7162	1 x 8ml			
Novolink™ DAB Substrate Buffer (Polymer)	RE7163	1 x 150ml			
Hematoxylin	RE7164	1 x 125ml			
Peroxidase Block	RE7165				1 x 5ml
Protein Block	RE7166				1 x 5ml
Post Primary	RE7167				1 x 5ml
Novolink™ Polymer	RE7168				1 x 5ml
DAB Chromogen	RE7169				1 x 1ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7171				1 x 5ml
Hematoxylin	RE7172				1 x 5ml

Reagentia	Productnr.	Novolink™ Max Polymer (1250 tests) RE7260-K	Novolink™ Polymer (250 tests) RE7200-K	Novolink™ Max DAB (Polymer) 1250 tests RE7270-K	Novolink™ DAB (Polymer) 250 tests RE7230-K
Post Primary	RE7111		1 x 25ml		
Novolink™ Polymer	RE7112		1 x 25ml		
Post Primary	RE7159	1 x 125ml			
Novolink™ Polymer	RE7161	1 x 125ml			
DAB Chromogen	RE7105				1 x 3ml
Novolink™ DAB Substrate Buffer (Polymer)	RE7143				1 x 30ml
DAB Chromogen	RE7162			1 x 8ml	
Novolink™ DAB Substrate Buffer (Polymer)	RE7163			1 x 150ml	

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Leica Biosystems Newcastle Ltd