

CytoVision* Version 7.4 - Intended Use

The CytoVision* Image Analysis and Capture System is a rapid metaphase finder, image acquisition and computer aided chromosome analysis system which assists the operator in viewing chromosomes and looking for cellular anomalies. CytoVision* enables a qualified Cytogeneticist to rapidly and accurately analyze the chromosome banding pattern.

All diagnostic decisions are made by the qualified clinician.

CytoVision* is approved for In-Vitro Diagnostic Use for CEP XY and Karyotyping.

CytoVision* CEP XY_ENG is an accessory to the CEP® X Spectrum Orange™/CEP® Y Spectrum Green™ DNA Probe kit (Abbott Laboratories, Abbott Park, Illinois, U.S.A) and is limited to the analysis of CEP XY probes via high magnification capture and analysis of interphase nuclei. CEP XY is indicated for use to assess the effectiveness of bone marrow transplantation in opposite-sex transplants.

Rx Only

CytoVision* is the trademark of Leica Biosystems Richmond, Inc. All third party trademarks are the property of their respective owners.

Leica Biosystems reserves the right to amend specifications without notice. All Rights Reserved. © Leica Biosystems Richmond, Inc. 2013



Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
U.S.A
Tel (800)-537-4669



CEpartner4U,
ESDOORNLAAN 13,
3951DB, MAARN,
THE NETHERLANDS,
www.cepartner4u.eu

*Reg. US Pat & TM Off and in other jurisdictions throughout the world