



Aperio VERSA Oil Scanning – Look Closer!

Enhance your Aperio VERSA 8 or 200 slide scanner with an Automated Oiler & High Resolution Oil Immersion Lens

Leica Biosystems is pleased to announce the release of the Aperio VERSA V2 Automated Oiler.

With this release, we are extending the capabilities of the Aperio VERSA system to include automated slide oiling and scanning with high resolution oil immersion objectives.


Aperio VERSA Oil Scanning Benefits

- ›› Lower exposures required for imaging
- ›› Higher transmission of signal
- ›› Greater image clarity
- ›› Improved intracellular detail & texture
- ›› 40x & 63x oil immersion objectives available

Start looking closer, with Aperio VERSA.
Contact your Leica Biosystems representative today!

For Research Use Only.
Not for use in diagnostic procedures.

LeicaBiosystems.com



BCL2 FISH at 63x oil



Pap smear at 40x oil