


Novocastra™ IHC Diluent

Product Code: RE7189

Intended Use	FOR RESEARCH USE ONLY.
Introduction	Novocastra™ IHC Diluent is intended for use as a diluent for Novocastra™ monoclonal and polyclonal primary antibodies. Do not use Novocastra™ IHC Diluent for the reconstitution of lyophilized antibodies.
Reagents Provided	Novocastra™ IHC Diluent - 500 mL
Storage and Stability	Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the product label.
Warnings and Precautions	 The concentration of ProClin™ 950 is 0.35%. It may cause irritation to the skin, eyes, mucus membranes and upper respiratory tract. A Material Safety Data Sheet is available upon request or available from www.LeicaBiosystems.com .
Procedure	Mix appropriate volumes of primary antibody and Novocastra™ IHC Diluent to achieve the optimal dilutions. Novocastra™ IHC Diluent can also stabilize diluted primary antibodies when stored at 2–8 °C.
Bibliography	Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

